
Characterizing and Developing the Clinical Grade Next Generation Sequencing based Gut Microbiome Assay with the Bioinformatics Solution


Abstract ID. 4704

Danyi Wang¹, Brajendra Kumar², Aaron Tenney³, Ravi Kiron⁴, Yang Liu³, Fei Zhong³, Juergen Scheuenpflug⁵, Zheng Feng ^{1#},

1. Clinical Measurement Sciences, Global Research & Development, EMD Serono, Billerica, MA, USA, A business of Merck KGaA, Darmstadt, Germany 2. Sigma-Aldrich Chemicals Pvt. Ltd., Bangalore, India (An affiliate of Merck KGaA, Darmstadt, Germany) 3. Bioinformatics R&D, Information Technology, Millipore Sigma, St. Louis, MO, USA 4. Strategic Innovation, Silicon Valley Innovation Hub, EMD Serono, San Jose, CA, USA 5. Merck KGaA, Darmstadt, Germany # Corresponding author




Taxonomic Profiles Varies among Extraction Methods




- The PCoA (Principal Coordinate Analysis) 3D visualization of the Weighted Unifrac beta diversity shows that kits A and C cluster closely together while kit B appears as an outlier. • The kit A sequencing samples cluster more closely together than both the other kits. • The taxonomic profiles of kit B have lower recall when compared to the known mixture profiles

- indicating that kit B was inefficient at detecting some of the bacteria.

RESULTS – Phase 1b




Dis(similarity) metrics-based ordination analysis for PCoA




Axis 3 (6.742 %) • Beta diversity distant metrices (Bray Curtis, Jaccard & Un-weighted Unifrac) based group significance test using pairwise PERMANOVA indicates no significant group separation between the kits. • P-value and the F-ratio value of pairwise permanova indicates the significant group separation based upon the based on weighted unifrac (P = 0.046)

Weighted Unifrac metric-based Beta group significance





pseudo- F ratio value in Pairwise Permanova results indicates that Kit C species abundance is significantly diverged from both Kit A and Kit B respectively.


Alpha diversity measures indicates richness and evenness within the kits



p-value and H-test does not show the statistically significant differences between groups







This study was sponsored by EMD Serono, a business of Merck KGaA, Darmstadt, Germany.



Relative frequency-based group-sample level composition

ACKNOWLEDGEMENTS AND DISCLOSURES